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Dear Customer,

Please read through these instructions for use carefully and if you have any questions, contact your physician or medical retailer or visit info@bauerfeind.com.

Intended use

Orthopedic orthoses are functional orthoses for the feet and are suitable for use in treating and preventing foot complaints. The ErgoPad ski & skate also makes ski boots more comfortable to wear and helps stabilize the movement of the foot during sporting activities. Suitable for a range of winter sports, particularly downhill and cross-country skiing, ice hockey and inline skating.

Indications

The ErgoPad ski & skate helps prevent and treat overloading complaints associated with specific types of sport. These include:

- Non-specific overloading complaints involving the foot
- Heel pain and talalgia
- Achilles tendon complaints
- Tibia complaints

It is also used to prevent and treat sports-related complaints experienced by patients with mild to severe arch decreased/splay foot or mild pes cavus/high arched foot.

Risks of using this product

⚠ Caution: Before using/applying a foot orthosis, you should consult a physician or specialist on the risks involved.

If the available internal volume of the shoe is not sufficient for the foot and the selected foot orthosis, or if the foot orthoses are worn in shoes of an incorrect length and width (or shoes with inadequate lateral guidance in the case of open shoe types), problems relating to poor fit may occur (e.g. a feeling of constriction or pressure, rubbing).

Please remember that shoes and foot orthoses form a single functional unit. If the feet are positioned on the edges of the foot orthosis, blisters and calluses may form on the soles of the feet. The same symptoms may occur if excessive mechanical strain leads to holes or creases developing in the top cover and/or cushioning material, or causes them to break. Foot orthoses that are damaged in this way must no longer be worn!

Contraindications

People with diabetes or rheumatoid foot deformities must seek medical advice before using this product to ensure optimal treatment.

Application instructions

Shoes suitable for wearing with foot orthoses must have sufficient internal volume for both the foot and the selected foot orthosis (adequate shoe length/width and, if necessary, a removable, sufficiently thick cushioned sole, as well as lateral guidance for the foot orthoses in the case of open shoe types). An uncomfortable fit that constrains the foot must always be avoided.

When selecting the correct size of foot orthoses, the packaging information serves merely as a guide. Only the correctly fitting shoes; foot orthoses must be at least 1.0 cm longer than the foot when the foot is in a weight-bearing, standing position. Please check the length allowance for both feet, as they may differ in length. In case of doubt, base your size selection on the longer foot.

The foot orthosis must be positioned straight in the shoe and have as much surface contact as possible with the shoe insole. Contact between the foot orthosis and the shoe upper helps reduce the risk of pressure marks or blistering on the skin of the soles of the feet.

Foot orthoses must always be worn as a pair! Any exceptions to this must be stipulated by the treating physician.

The ErgoPad ski & skate is designed to be worn repeatedly by one user only.

Warm-up: Foot orthoses are items for daily use and as such they are subject to natural wear when used on a regular basis. Slight discoloration, in particular fading of the top covers of the foot orthoses, may occur during their use. This is dependent on a number of factors, such as moisture levels and sock/stocking materials. This is not covered by warranty.

Warranty is provided in accordance with the provisions as well as the points mentioned and contractual terms of Bauerfeind AG (Terms and Conditions), as long as all the terms stipulated in these instructions for use are observed.

Disposal

You can dispose of the product in accordance with national legal requirements.

Modification instructions

Fitting: For an optimum fit, there should be a length allowance of approximately 1 cm between the foot and the foot orthosis and / or shoe, in a standing position. If the foot orthoses are not the right size for your shoes, you can cut them down or reshape them using thermforming. Please follow the instructions below.

Cutting: Take the inlays out of your downhill skiing / cross-country skiing / ice hockey boots (or the inner shoes where applicable) and place them on your ErgoPad ski & skate.

Align the heel sections so that they are positioned precisely on top of one another and draw around the inlay on your ErgoPad ski & skate.

Cut along the outline you have drawn in the forefoot area. (Make sure you do not cut the foot orthosis too small. We recommend cutting off less material than you need to at first and trying on the foot orthosis before carrying on cutting).

Place the ErgoPad ski & skate in your shoes (foot orthoses for left/right feet are marked as such on the back).

Thermoforming

In exceptional cases and having sought specialist advice, it may be necessary to use thermoforming to adapt the foot orthoses to your foot.

1. After cutting down your ErgoPad ski & skate, preheat your oven to 90 °C/200 °F.
2. Put the foot orthoses (with the top covers facing forward on a clean surface) in the oven for approximately 2 minutes.


Risk of burning! Wear protective gloves when taking the foot orthoses out of the oven.

As soon as you have taken the foot orthoses out of the oven, place them immediately in your shoes. Do not put on your shoes straight away (make sure you are wearing socks or stockings).

5. Then stand for 2 minutes with your feet shoulder-width apart and pointing forward. Your foot orthoses should now be adapted to the shape of your feet.

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Material	Material
Polyester (PES)	Polyester (PES)
Polyurethan (PUR)	Polyurethane (PUR)

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Bauerfeind Nordic AB