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en english

## Dear Customer,

Please read through these instructions for use carefully and if you have any questions contact your physician, medical retailer or our Customer Service department (info@bauerfeind.com).

### Uses / Indications

This foot orthosis is suitable for the conservative treatment of orthopedic foot and knee complaints and some hip and back complaints.

Orthopedic foot orthoses are used in the treatment of foot complaints and foot deformations, such as flat foot, splay foot, pes valgus, mild cases of pes cavus, or a combination of these conditions. The treating physician is responsible for the precise diagnosis.

Orthopedic orthotic blanks must be properly fitted by a qualified specialist<sup>1</sup>.

In addition to the treatment of foot complaints, foot orthoses can also be used to prevent excessive strain on the entire musculoskeletal system (knee, hip and back regions).

### Risks of using this product

⚠ If the available internal volume of the shoe is not sufficient for the foot and the selected orthosis, or if the foot orthoses are worn in shoes which are the incorrect length and width (or shoes with inadequate lateral guidance in the case of open shoe types), problems relating to poor fit may occur (e.g. a feeling of constriction or pressure, rubbing).

If the feet are positioned on the edges of the foot orthosis, blisters and calluses may form on the soles of the feet. The same symptoms may occur if excessive mechanical strain leads to holes or creases developing in the top covers, cushioning material or orthotic cores, or to the breakage of these. Foot orthoses which are damaged in this way must no longer be worn!

Creases in stockings or pantyhose must be pulled flat before the shoes are put on, as creases of this kind can cause blisters and pressure points.

If your symptoms increase, please contact your physician or medical retailer without delay.

The foot orthosis may not be used to treat patients with a diabetic foot syndrome or rheumatoid foot deformities.

### Application instructions

The ErgoPad work for safety shoes must be fitted to the feet and / or shoes by a qualified specialist.

When selecting the correct size of foot orthosis, the packaging information serves merely as a guide. Just like correctly fitting shoes, orthoses must be at least 1.0 cm longer than the foot, when the foot is in a weight-bearing, standing position **1**.

The shoes must guarantee sufficient internal volume for both the foot and the selected foot orthosis (adequate shoe length / width, removable, sufficiently thick cushioned sole, adequate lateral guidance of the foot orthoses in the case of open shoe types). An uncomfortable fit that constricts the foot must always be avoided **2**.

The forefoot (a1), midfoot (a2) and hindfoot (a3) must be supported sufficiently evenly from below by the foot orthosis and must not be positioned on one of its edges **3**.

Foot orthoses must always be worn as a pair! Any exceptions to this must be stipulated by the treating physician.

It must always be ensured that the left and right foot orthoses are clearly identified in a way which is easily understood by the patient.

The foot orthoses must not be exposed to temperatures exceeding 60 °C (e.g. heating, hair-dryer, direct sunlight).

The materials used for the foot orthoses are subject to aging and their properties can therefore change.

### Storage conditions

The products must not be exposed to direct UV light and must be stored within a temperature range of 10 °C to 25 °C and free from any mechanical influence. The

**USA**  
Bauerfeind Ges.m.b.H. Hainburger Straße 33 1030 Wien P +43 (0) 800 4430-130 F +43 (0) 800 4430-130 E info@bauerfeind.at

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manufacturer recommends that regular checks are carried out on both the foot orthoses in storage and the storage conditions.

### Instructions for cleaning and disinfection

Foot orthoses must be removed from the shoes on a daily basis in order to allow moisture to evaporate.

All surface materials facing the foot can be wiped if required, using a mild cleanser and a sponge or cloth. After cleaning, any residual soap must be thoroughly rinsed off with water. The foot orthoses cannot be washed in the washing machine.

The orthoses must be left to air dry and no external heat source should be used (heating, hair-dryer, direct sunlight).

Mild disinfectants suitable for the human skin may be applied sparingly (sprayed on). However, in this case, discoloration of the top cover materials may occur.

### Permissible conditions and places of use

The foot orthoses must be worn in conjunction with hosiery in suitable shoes with a sufficient internal volume for both the foot and the selected orthosis.

### Notes on usage

The product is designed to provide care for an individual patient / individual orthosis wearer and must not be reused or used by any other person. The manufacturer accepts no liability in the event of improper use of the product.

Foot orthoses are affected by various mechanical forces (body weight, tension, pressure, torsion, friction), with the result that their serviceable life depends on the load. Temperatures, stocking materials, the properties of the orthosis materials, shoes and environmental influences (e.g. perspiration) all play their part in reducing the functional life of the foot orthosis.

After a foot orthosis has been worn for a maximum of six months, the manufacturer recommends that the orthopedic functions of the orthosis are checked by a qualified specialist.

### Warranty

The statutory regulations of the country of purchase apply. Where relevant, the country-specific warranty regulations between the buyer and retailer apply. Please contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty.

Please do not make any alterations to the product yourself. This particularly applies to products that have been individually adjusted by a trained specialist. Please follow our instructions for use and care.

They are based on years of experience and ensure that our medically effective products continue to work for a long time. After all, only aids that function perfectly will help you. Failure to follow the instructions can also limit the warranty.

### Duty to report

Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of a medical product to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

### Disposal

You can dispose of the product in accordance with national legal requirements.

### Modification instructions for fitting the ErgoPad work for safety shoes in accordance with DIN EN ISO 20345 and DIN EN 61340-5-1

The following materials must be used for further modification and bonding purposes:

<b>Pryx® antistatic</b> Purchase details: Fritz Minke GmbH & Co KG Cecilienstraße 31 D-47051 Duisburg Phone +49 (0) 203 28101-0 Fax +49 (0) 203 28101-70 E-mail gminke@minke.de	<b>Renol =AL= (antistatic)</b> Item No. 24002601 (640g) Purchase details: Ortho Planet Leder Brinkmann GmbH Eduard-Pestel-Str. 3 D-49080 Osnabrück Phone +49 (0) 541 95933-0 Fax +49 (0) 541 95933-33 E-mail info@leder-brinkmann.de
--	---

### Mounting and modification specifications

The foot orthosis may only be modified using the procedure

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