



Certificate

for the company

Bauerfeind AG
Triebeser Straße 16
07937 Zeulenroda-Triebes

manufactured products groups
(listed in the annex)

can be certified after testing
Cytotoxicity with the neutral red method and the method of LDH-release
according to EN ISO 10993-5:2009 and
Irritation according to EN ISO 10993-10:2013,

the biological harmlessness according to EN ISO 10993-1:2020.

HygCen Germany GmbH | Bornhövedstrasse 78 | 19055 Schwerin

Identification number:

99430

Schwerin, 2022-08-31

Valid until:

2025-08-31

DocuSigned by:

A handwritten signature in blue ink, appearing to be 'S. Werner', is enclosed within a blue rounded rectangular border.

BD051B9D54C746F...

Dr. med. univ. Sebastian Werner
Head of Scientific-Technical Affairs

The following test samples were tested according to EN ISO 10993-1:

Globofit (PUR ca. 12° Shore A):

beige (previously haut) - 104961, 104962 - SN 14938
hellblau - 106719, 104960, 106961 - SN 9378
Plastazote beige (previously haut) - 104953 - SN 9392

Lunairflex (EVA 22° ShoreA):

beige (previously haut) - 105006, 105007, 105164, 109402, 108410 - SN 34310

Lunasoft SL (EVA 40° Shore A):

mehrfarbig - 105014, 105015 - SN 9381

Poron combination (PUR, EVA/PE):

lachsfarben - 106960, 107202 - SN 30355

TECfom (EVA/PE 25° Shore A):

blau und saphirblau - 107306, 107870, 107880, 107881 - SN 9383
gold - 113214, 115692, 115711 - SN 24373, SN 34833, SN 34896
grau - 107882-107884 - SN 22878
weiß - 107918 - SN 34971
beige (previously haut) - 107919 - SN 34972
hellazurblau - 107037, 107872, 108000 - SN 9383
schwarz - 108062, 108063, 109677 - SN 34973

TEClast (EVA/PE 45° Shore):

blau - 110402 - SN 15168
grau - 108441, 110401, 112495, 113051 - SN 14931, SN 34300

TECmoll (EVA/PE 35° Shore A):

graubraun - 108061, 108064 - SN 23043
blau und saphirblau - 107857, 107902, 107873, 107886 - SN 9384
gelbgrün - 112502, 113990 - SN 22877
gold - 110511 - SN 22878
grau und silbergrau - 107853, 107874, 107855, 107792, 108350, 107792, 107801,
107874, 107889, 107890 - SN 34313

The following test samples were tested according to EN ISO 10993-1:

TECsoft (PE 25° Shore A):

grau - 109460-109461 - SN 13235

TECstrong (EVA/PE 50°-55° Shore A):

braun - 105001 - SN 12538

blau - 105003 - SN 9334

TECtop (EVA/PE 30° Shore A):

weiß - 104997, 107918 - SN 9386, SN 34971

arktisblau - 114516 - SN 34304

grau und silbergrau - 107875, 107910-107912, 107923, 107925 - SN 9397

orange - 114514 - SN 34303

rot - 104985 - SN 34301

schwarz - 107475, 108170, 109840 - SN 34312

terrabraun - 107927 - SN 9397

gelbgrün - 114515 - SN 34305

Other:

Nessel Roh Baumwolle - 105108 - SN 34309

Sporteva orange - 115359 - SN 33088

XVA snow white - 115780 - SN 35054